



PRODUCT CODE	INNER DIAMETER [μm] ± 2 %	OUTER DIAMETER [μm] ± 5 %	COATING DIAMETER [μm] ± 5 %	STANDARD COATING	SPECIFIC FEATURE
<b>POLYIMIDE COATED CAPILLARIES (-190°C TO 385°C)</b>					
CAP50/120PI	50	120	144	POLYIMIDE	
CAP50/330PI	50	330	363	POLYIMIDE	
CAP75/120PI	75	120	144	POLYIMIDE	
CAP75/330PI	75	330	363	POLYIMIDE	
CAP100/120PI	100	120	144	POLYIMIDE	
CAP100/330PI	100	330	363	POLYIMIDE	
CAP127/160PI	127	160	190	POLYIMIDE	
CAP127/210PI	127	210	240	POLYIMIDE	
CAP150/330PI	150	330	363	POLYIMIDE	
CAP180/330PI	180	330	363	POLYIMIDE	
CAP200/270PI	200	270	300	POLYIMIDE	
CAP250/340PI	250	340	363	POLYIMIDE	
CAP300/400PI	300	400	430	POLYIMIDE	
CAP350/470PI	350	470	500	POLYIMIDE	
CAP400/540PI	400	540	570	POLYIMIDE	
<b>UNCOATED CAPILLARIES (TO 1000°C)</b>					
CAP128/500	128	500			
CAP128/900	128	900			
CAP550/700	550	700			
CAP550/900	550	900			
CAP800/1000	800	1000			
CAP800/1320	800	1320			
CAP950/1200	950	1200			
CAP950/1570	950	1570			

FOR UV-VIS OPTICAL QUALITY, ADD **UV** TO THE PRODUCT CODE.

FOR VIS-NIR OPTICAL QUALITY, ADD **IR** TO THE PRODUCT CODE.

**OTHER SPECIFICATIONS UPON REQUEST.**